

# COMPASS LE OWC Root Exclusion: Porewater Nutrients

Apr-May 2024 Samples

2025-10-30

## Contents

0.1	Run Information . . . . .	2
0.2	Plot standard curves . . . . .	3
0.3	Check NOx Reduction Efficiency . . . . .	9
0.4	Analyze the Check Standards . . . . .	10
0.5	Analyze Blanks . . . . .	11
0.6	Check to see if samples run match metadata & merge info . . . . .	12
0.7	Visualize data . . . . .	13

## 0.1 Run Information

```
#set the run date & user name
run_date <- "08/20/2025"
sample_year <- "2024"
sample_month <- "Apr-May"
user <- "Roberta Peixoto"
QAQC_by <- "Zoe Read"

#identify the files you want to read in
#read in as a list to accommodate multiple runs in a month
data_files <- c("Raw Data/Seal_COMPASS_LE_OWC_RootExclusion_NUTR_Apr-May2024.csv")

# Define the file path for QAQC log file - NO Need to change just check year
file_path <- "Raw Data/SEAL_COMPASS_LE_OWC_RootExclusion_QAQC_Log_2025.csv"
final_path <- "Processed Data/COMPASS_LE_OWC_RootExclusion_Nutrients_2024_Apr-May2024.csv"

#record any notes about the run or anything other info here:
run_notes <- "The last CCV check std for Phosphate is very high.
All other check standards are too low. " #Add notes from QAQC

#Set up file path for metadata
#downloaded metadata csv - downloaded from Google drive as csv for this year
Raw_Metadata = "Raw Data/COMPASS_OWC_Root_Exclusion_Nutrient_SampleLog.csv"

cat(run_notes)
```

```
## The last CCV check std for Phosphate is very high.
```

```
## All other check standards are too low.
```


```
##Setup
```

```
##Read in metadata and create similar sample IDs for matching to samples
```

```
##Import Data & Clean
##Assessing standard Curves
```


## 0.2 Plot standard curves

```
## 'geom_smooth()' using formula = 'y ~ x'
```




```
## 'geom_smooth()' using formula = 'y ~ x'
```

## NH3 Standard Curve




```
## 'geom_smooth()' using formula = 'y ~ x'
```

# PO4 Standard Curve




## SO4 Standard Curve



```
## 'geom_smooth()' using formula = 'y ~ x'
```

## NO2 Standard Curve



```
## [1] "NOx Curve r2 GOOD - PROCEED"

## [1] "NH3 Curve r2 GOOD - PROCEED"

## [1] "PO4 Curve r2 GOOD - PROCEED"

## [1] "SO4 Curve r2 GOOD - PROCEED"

## [1] "NO2 Curve r2 GOOD - PROCEED"

## [1] "QAQC log file exists and has been read into the code."

## 'geom_line()': Each group consists of only one observation.
## i Do you need to adjust the group aesthetic?
```




Table 1: Average Slope by Analyte


Test	avg_slope
NH3	18.657
NO2	437.467
NOx	1.252
PO4	2.616
SO4	1.570



##Dilution Corrections - ensure the latest dilution is kept

## No duplicated samples.

### 0.3 Check NO<sub>x</sub> Reduction Efficiency



## [1] "Mean NO<sub>x</sub> Reduction Efficiency >95% - PROCEED"

## [1] 100.045

## 0.4 Analyze the Check Standards


```
## [1] "NOx Check Standard RSD within Range - PROCEED"
```

```
## [1] "NH3 Check Standard RSD within Range - PROCEED"
```

```
## [1] "PO4 CHECK STANDARD RSD TOO HIGH - REASSESS"
```

```
## [1] "SO4 Check Standard RSD within Range - PROCEED"
```

```
## [1] "NO2 Check Standard RSD within Range - PROCEED"
```



```
## [1] "<60% of NOx Check Standards are within range of expected concentration - REASSESS"
```

```
## [1] "<60% of NH3 Check Standards are within range of expected concentration - REASSESS"
```

```
## [1] "<60% of PO4 Check Standards are within range of expected concentration - REASSESS"
```

```
## [1] "<60% of SO4 Check Standards are within range of expected concentration - REASSESS"
```

```
## [1] "<60% of NO2 Check Standards are within range of expected concentration - REASSESS"
```

## 0.5 Analyze Blanks

## [1] ">60% of NOx Blanks are below the lower 25% quartile of samples or 1/2 detection limit - PROCEED"

## [1] ">60% of NH3 Blanks are below the lower 25% quartile of samples - PROCEED"

## [1] ">60% of PO4 Blanks are below the lower 25% quartile of samples- PROCEED"

## [1] ">60% of SO4 Blanks are below the lower 25% quartile of samples- PROCEED"

## [1] ">60% of NO2 Blanks are below the lower 25% quartile of samples- PROCEED"




Table 2: Mean Concentration of Blanks

Test	Blank_Mean_Conc
NOx	-0.0029
NH3	-0.0037
PO4	0.0022
SO4	0.0433
NO2	-0.5667


```
##Sample Flagging - Within range of standard curve  
##Pull out sample id information
```

## 0.6 Check to see if samples run match metadata & merge info


```
## Check Sample IDs with Metadata
```

```
## All sample IDs are present in metadata.
```


```
## Visualize Data
```





## Porewater NH3



## Porewater PO4





##Export Processed Data

#end